

## Intended Use:

This respirator is intended to protect the user of certain airborne particles. Provides respiratory protection to the wearer greater than 95%. This product is disposable and it can fit a wide range of face sizes.

## Important:

Before using the respirator, follow the instructions for use and adjustment.

## Description:

Respirator N95 Particulate filter (>95% filter efficiency level) effective against particulate aerosols free of oil; time use restrictions may apply.

## Instructions for Use

Follow these instructions each time the respirator is worn. (See Figures 1-5)

- 1 Hold the respirator in your hand with the elastic straps hanging down.
- 2 Place the respirator over your nose and mouth.
- 3 Place the top strap over your head and place the bottom elastic strap behind your neck.
- 4 Carefully adapt the nasal clip to the shape of the nose so that it fits comfortably and adjust the respirator to achieve a completely tight fit around the face. Check the tightness of the facial fit of the respirator with a puff.
- 5 If air leaks occur in the nasal area (nasal clip) or elsewhere you must remove them by readjusting the nasal clip and respirator the contaminated area or risk zone.



## Cautions and Limitations

- Not for use in atmospheres with less 19% oxygen.
- Not for use in atmospheres that represent a danger to life or health.
- Should not be used in places with higher concentrations of pollutants than those established by regulated standards.
- Incorrect use and maintenance of this product could result in illness or death.
- All approved respirators shall be selected fitted used and maintained according to the applicable country regulations.
- You must not replace, modify, or add parts.
- Refer to the instructions for use and / or maintenance manuals for proper use of the respirator.
- Before using the product, please refer to the instructions for the user.

## Time Use Limitation


If breathing becomes difficult or the respirator becomes damaged or soiled, leave the contaminated area immediately and replace the respirator.

## Storage Conditions and Shelf Life

Shelf life of unopened product is 5 years from date of manufacture when stored within temperature range of -20°C to +30°C and at less than 80% relative humidity.



Shelf life of unopened product is 5 years from date of manufacture.



Temperature Range: -20°C To +30°C



Maximum Relative Humidity: <80%

